

# A Novel Multiplex PCR Assay with Melting Curve Analysis for Rapid Detection of *Aspergillus* Species and Azole Resistance

Shu Li Li<sup>1,2</sup>, Zhi Xian Wang<sup>1,2</sup>, Doudou Jiang<sup>1,2</sup>, Ying Wang<sup>1,2</sup>, Yuan Zhang<sup>1,2</sup>, He Wang<sup>1,2</sup>, Yan Su<sup>1,2</sup>, Ze Qi Zhou<sup>1,2</sup>

<sup>1</sup>Dynamiker Sub-Center of Beijing Key Laboratory for Mechanisms Research and Precision Diagnosis of Invasive Fungal Disease; <sup>2</sup>Tianjin Enterprise Key Laboratory for precision diagnosis technology of invasive fungal diseases

## BACKGROUND

*Aspergillus* spp. are filamentous, environmental fungi that cause a wide spectrum of infections in humans, including hypersensitivity reactions, chronic pulmonary infections, and acute life-threatening infections. The diagnosis of aspergillosis remains difficult, and the rapid emergence of azole resistance in *A. fumigatus* is worrisome. The novel multiplex PCR assay with melting curve analysis provides a solution in this matter. The aim of this study was to validate the MycoMDx *Aspergillus* PCR Assay Kit and *A. fumigatus* resistance multiplex assay Kit produced by Dynamiker Biotechnology (Tianjin) Co., Ltd., and to evaluate its performance on clinical samples for the detection of four common *Aspergillus* species (*A. fumigatus*/*A. flavus*/*A. terreus*/*A. niger*) and three resistance-associated mutations (L98H/Y121F/T289A) in the Cyp51A gene. (Figure 1).

## METHOD

The diagnostic criteria of this study were based on the 2019 EORTC/MSGERC Definitions of invasive fungal diseases. Thirty-two patients with proven/probable aspergillosis and fifty non-aspergillosis patients were enrolled and detected for *Aspergillus* PCR and azole resistance.

## RESULT

The results showed that the sensitivity and specificity of *Aspergillus* PCR in bronchoalveolar lavage fluid (BALF) were 93.75% (30/32) and 98.00% (49/50). A total of four common clinical *Aspergillus* species were detected, among which *A. fumigatus* appeared the most frequently, followed by *A. flavus*, *A. terreus*, *A. niger*. In addition, resistance-associated mutations were detected in 6 BALF samples by *A. fumigatus* resistance multiplex assay, including 3 samples containing L98H mutation, 2 samples containing Y121F mutation, 1 sample containing T289A mutation.

## CONCLUSION

The results of this study indicate that the novel multiplex PCR assay with melting curve analysis can identify four medically important *Aspergillus* species DNA and detect presence of the L98H/Y121F/T289A mutation.

## KEY WORDS

Multiplex PCR; Azole Resistance; *Aspergillus* species




Figure 1 The workflow of a novel multiplex PCR assay