

Diagnostic Value of Combined Detection of BDG assay and GM assay in Invasive Aspergillosis

Poster Number: A-328

Mei Li¹, Shigui Liu¹, Wang Yao¹, Yuan Zhang¹, He Wang¹, Yan Su¹, Zeqi Zhou¹

¹ Tianjin Enterprise Key Laboratory for precision diagnosis technology of invasive fungal diseases - Tianjin (China)



丹娜生物
DMK



INTRODUCTION

- With the increasing number of patients who are immunosuppressed or immunocompromised, there has been a corresponding rise in invasive fungal infections over the past few decades, which are associated with high morbidity and mortality rates, ranging from 60% to 90%. Therefore, the role of early and accurate diagnosis is crucial.
- Conventional fungal identification relies on microscopic analysis and the growth of microorganisms in culture media. In contrast, serological testing offers a faster method for detecting the causative fungi, thereby aiding in the diagnostic decision-making process. The advantages of serology-based tests include the rapid results obtained, in contrast to culture methods, and the non-invasive nature of the samples (e.g., blood, urine, sputum), while also serving as potential prognostic markers, such as (1, 3)- β -D-glucan (BDG assay) and galactomannan (GM assay).
- Recently, a novel chemiluminescent immunoassay (CLIA) for BDG and GM developed by Dynamiker Biotechnology (Tianjin) Co., Ltd. was introduced as a screening test for invasive aspergillosis (IA, **Figure 1**). This fully automated assay offers high throughput and a short detection time.

AIM

In this study, we explored the diagnostic value of the combined detection of BDG-CLIA and GM-CLIA for IA.

CONCLUSIONS

- We conclude that BDG-CLIA demonstrates modest performance in identifying IA patients.
- GM-CLIA exhibits superior accuracy compared to BDG-CLIA.
- **The combination of both improves the sensitivity at the cost of specificity, and can be used as an auxiliary diagnosis for IA.**




Figure 1 Detection Principle of the Dynamiker Fungus (1,3)- β -D Glucan Assay (Chemiluminescence, A) and Dynamiker Aspergillus galactomannan Assay (Chemiluminescence, B); Fully Automated Chemiluminescence Analyzer (C).

METHOD

A total of 114 patients were admitted and enrolled from March 2023 to July 2023. According to the 2019 EORTC/MSG definitions of invasive aspergillosis, the patients were divided into two groups: the IA group (67 cases) and the non-IA group (47 cases). **The samples were tested in parallel using the Dynamiker BDG-CLIA and GM-CLIA assays.**

RESULTS

Compared to using BDG-CLIA alone (sensitivity 68.66%) or GM-CLIA alone (sensitivity 86.57%) the combination strategy demonstrated a significantly higher sensitivity (92.53%, $p < 0.01$) in the IA group (Table 1).

Table 1 Diagnostic value of combined detection of BDG-CLIA and GM-CLIA in IA

	Sensitivity [% (95% CI)]	Specificity [% (95% CI)]	Positive predictive value [% (95% CI)]	Negative predictive value [% (95% CI)]
BDG-CLIA	68.66 (56.03-79.13)	70.21 (54.92-82.21)	76.67 (63.66-86.22)	61.11 (46.87-73.76)
GM-CLIA	86.57 (75.53-93.30)	82.98 (68.65-91.86)	87.88 (76.96-94.25)	81.25 (66.90-90.56)
Combined	92.53 (82.74-97.22)	68.09 (52.75-80.48)	80.52 (69.60-88.33)	86.49 (70.43-94.92)