



Doudou Jiang^{1,2}, Shuli Li^{1,2}, Zhixian Wang^{1,2}, Ying Wang^{1,2}, Yuan Zhang^{1,2}, He Wang^{1,2}, Yan Su^{1,2}, Zeqi Zhou^{1,2}

¹Dynamiker Sub-Center of Beijing Key Laboratory for Mechanisms Research and Precision Diagnosis of Invasive Fungal Disease - Tianjin (China)

²Tianjin Enterprise Key Laboratory for precision diagnosis technology of invasive fungal diseases - Tianjin (China)



Background

- Under the selective challenge of antibiotics, the variety and number of drug-resistant pathogenic microorganisms have increased significantly, which brings great challenges to clinical diagnosis and treatment, especially the infection caused by carbapenem resistant *Enterobacteriaceae* (CRE).
- The production of carbapenemases, which are members of the β -lactamases that are divided into classes A, B and D based on their molecular structures, is the major mechanism underlying carbapenem resistance around the world and represents a great health concern. These are *Klebsiella pneumoniae* enzyme (KPC), New Delhi metal β -lactamase (NDM), carbapenem hydrolyzed oxalase (OXA-48 like), integrin-encoded metal β -lactamase (VIM) and IMP (para imipenem).
- The real-time fluorescence quantitative PCR (qPCR) method based on molecular beacons was combined with the melting curve analysis to identify five drug resistance genes simultaneously by a single PCR reaction, with rapid detection, high sensitivity and strong specificity.
- We developed the novel fluorogenic assay for rapid detection of Carbapenemases in multidrug-resistant *Enterobacteriaceae* (Dynamiker Biotechnology (Tianjin) Co., Ltd.).



Methods

- We evaluated the performance of the novel fluorogenic assay for rapid detection of Carbapenemases in multidrug-resistant *Enterobacteriaceae*, including the limit of detection (LoD) and cross-reactivity, and compared it with the lateral flow immunochromatography assay (LFA).



Results

- The LoD ranged from 75-450 CFU/mL for the five carbapenemase genes.
- Comparison of qPCR and LFA results from twenty-three CRE clinical isolates with characterized carbapenemase content demonstrated a complete agreement.
- The analytical specificity for target genes was 100%, as assessed with a panel of 15 pathogens, which indicated no cross-reactions.

Table 1 Detection results of clinical isolates by the two methods between qPCR and LFA

Strain No.	Organism name	PCR	LFA	Consistency
CRE 1-1	<i>Klebsiella pneumoniae</i>	VIM-1	VIM	Y
CRE 1-2	<i>Pseudomonas aeruginosa</i>	VIM-10	VIM	Y
CRE 2-1	<i>Escherichia coli</i>	NDM-5	NDM	Y
CRE 2-2	<i>Klebsiella pneumoniae</i>	NDM-1	NDM	Y
CRE 3-1	/	KPC-2/NDM-1	KPC/NDM	Y
CRE 3-2	/	KPC-2/NDM-1	KPC/NDM	Y
CRE 4-1	<i>Enterobacter cloacae</i>	IMP-4	IMP	Y
CRE 4-2	<i>Klebsiella pneumoniae</i>	IMP-4	IMP	Y
CRE 5-1	<i>Klebsiella pneumoniae</i>	OXA-48	OXA-48	Y
CRE 5-2	<i>Klebsiella pneumoniae</i>	OXA-232	OXA-48	Y
CRE 6-1	<i>Klebsiella pneumoniae</i>	KPC-2	KPC	Y
CRE 6-2	<i>Serratia marcescens</i>	KPC-2	KPC	Y
CRE-0	<i>Serratia marcescens</i>	/	/	Y
CRE-100	<i>Klebsiella pneumoniae</i>	/	/	Y
CRE-102	<i>Klebsiella pneumoniae</i>	KPC	KPC	Y
CRE-121	<i>Klebsiella pneumoniae</i>	NDM	NDM	Y
CRE-134	<i>Escherichia coli</i>	/	/	Y
CRE-144	<i>Klebsiella pneumoniae</i>	/	/	Y
CRE-152	<i>Escherichia coli</i>	NDM	NDM	Y
CRE-176	<i>Klebsiella pneumoniae</i>	NDM	NDM	Y
CRE-227	/	NDM	NDM	Y
CRE-231	/	VIM	VIM	Y
KM174	/	NDM	NDM	Y

"Y": Yes

Table 2 The cross-reaction results of 15 pathogens

Organism name	Strain No.	Result
<i>Escherichia coli</i>	ATCC8739	N
<i>Achromobacter baumannii</i>	ATCC17912	N
<i>klebsiella pneumoniae</i>	ATCC700603	N
<i>Enterobacter cloacae</i>	ATCC13047	N
<i>Candida albicans</i>	ATCC10231	N
<i>Klebsiella oxytoca</i>	ATCC15328	N
<i>Enterococcus faecalis</i>	ATCC29212	N
<i>Staphylococcus aureus</i>	ATCC25923	N
<i>Clostridium difficile</i>	ATCC43255	N
<i>Bifidobacterium adolescentis</i>	ATCC15703	N
<i>Enterobacter aerogenes</i>	ATCC13048	N
<i>Helicobacter pylori</i>	ATCC43504	N
<i>Campylobacter jejuni</i>	ATCC33560	N
<i>Cryptococcus neoformans</i>	ATCC 32045	N
<i>Cryptococcus gattii</i>	ATCC MYA-4560	N

"N": Negative

Conclusions

- The novel fluorogenic assay for rapid detection of Carbapenemases in multidrug-resistant *Enterobacteriaceae* is an accurate and rapid method to identify KPC, NDM, VIM, IMP and OXA-48-like carbapenemases in the clinical microbiology laboratory, which can guide infection control programs to limit the spread of these organisms.